

Product Datasheet

Recombinant Human FGF10 Protein (orb1551033)

Description	Recombinant Human FGF10 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	FGF10
Note	For research use only
Purity	> 98 % by SDS-PAGE.
Source	KMP1454, Recombinant Human FGF10 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gln38-Ser208) of human FGF10 (Accession #NP_004456.1).
Expiration Date	6 months from date of receipt.