

Product Datasheet

Recombinant Human TNFRSF1A Protein, His Tag (orb1551039)

Description	Recombinant Human TNFRSF1A Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	TNFRSF1A
Note	For research use only
Purity	> 92 % by SDS-PAGE.
Source	KMP1448, Recombinant Human TNFRSF1A Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Thr211) of human TNFRSF1A (Accession #NP_001056.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.