

Product Datasheet

Recombinant Human C5a Protein (orb1551041)

Description	Recombinant Human C5a Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	C5a
Note	For research use only
Purity	> 94 % by SDS-PAGE.
Source	KMP1446, Recombinant Human C5a Protein is produced by E. coli expression system. The target protein is expressed with sequence (Leu679-Arg751) of human C5a (Accession #NP_001726.2).
Expiration Date	6 months from date of receipt.