

Product Datasheet

Recombinant Human IL13RA1 Protein, Fc & His Tag (orb1551069)

Description	Recombinant Human IL13RA1 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL13RA1
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1418, Recombinant Human IL13RA1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met 1-Thr 343) of human IL13RA1 (Accession #NP_001551.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.