

Product Datasheet

Recombinant Human ANGPTL3/Angiopoietin-like 3 Protein, His Tag (orb1551072)

Description	Recombinant Human ANGPTL3/Angiopoietin-like 3 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ANGPTL3
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1415, Recombinant Human ANGPTL3/Angiopoietin-like 3 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser17-Pro220) of human ANGPTL3/Angiopoietin-like 3 (Accession #NP_055310.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.