

## Product Datasheet

# Recombinant Human CD155/PVR Protein, Fc & His Tag (orb1551074)

<b>Description</b>	Recombinant Human CD155/PVR Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	CD155
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1413, Recombinant Human CD155/PVR Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Asp28-Asn343) of human CD155/PVR (Accession #NP_006496.4) fused with an Fc, 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.