

Product Datasheet

Recombinant Human IL-29 Protein, His Tag (orb1551086)

Description	Recombinant Human IL-29 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	IL29
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1401, Recombinant Human IL-29 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gly20- Thr200) of human IL-29 (Accession #NP_742152.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.