

Product Datasheet

Recombinant Human IL7R alpha Protein, Fc & His Tag (orb1551087)

Description	Recombinant Human IL7R alpha Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL7R alpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1400, Recombinant Human IL7R alpha Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Glu21-Gly236) of human IL7R alpha (Accession #NP_002176.2) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.