

Product Datasheet

Recombinant Human Frizzled-1 Protein, Fc & His Tag (orb1551089)

Description	Recombinant Human Frizzled-1 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Frizzled1
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1398, Recombinant Human Frizzled-1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gln73-His Tag253) of human Frizzled-1 (Accession #NP_003496.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.