

Product Datasheet

Recombinant human VEGF 165 Protein, His Tag (orb1551104)

Description	Recombinant human VEGF 165 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	VEGF 165
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1383, Recombinant human VEGF 165 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala27-Arg191) of human VEGF 165 (Accession #NP_001165097.1) fused with a 6xHis Tag at the N-terminus.
Expiration Date	6 months from date of receipt.