

Product Datasheet

Recombinant human Follistatin Protein, His Tag (orb1551106)

Description	Recombinant human Follistatin Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Follistatin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1381, Recombinant human Follistatin Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gly30-Trp344) of human Follistatin (Accession #NP_037541.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.