

Product Datasheet

Recombinant human Beta-catenin Protein, His Tag (orb1551109)

Description	Recombinant human Beta-catenin Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Betacatenin
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1378, Recombinant human Beta-catenin Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Leu781) of human Beta-catenin (Accession #XP_016861227.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.