

## Product Datasheet

# Recombinant human Beta-catenin Protein, His Tag (orb1551111)

<b>Description</b>	Recombinant human Beta-catenin Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	Betacatenin
<b>Note</b>	For research use only
<b>Purity</b>	> 90% by SDS-PAGE.
<b>Source</b>	KMP1376, Recombinant human Beta-catenin Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Leu781) of human Beta-catenin (Accession #XP_016861227.1) fused with a 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.