

Product Datasheet

Recombinant human GTPase KRas/KRAS Protein, His Tag (orb1551115)

Description	Recombinant human GTPase KRas/KRAS Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	GTPase KRas
Note	For research use only
Purity	> 85% by SDS-PAGE.
Source	KMP1372, Recombinant human GTPase KRas/KRAS Protein is produced by E. coli expression system. The target protein is expressed with sequence (Thr2-Cys186) of human GTPase KRas/KRAS (Accession #NP_203524.1) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.