

Product Datasheet

Recombinant human CD24 Protein, Fc His Tag (orb1551116)

Description	Recombinant human CD24 Protein, Fc His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD24
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1371, Recombinant human CD24 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Gly59) of human CD24 (Accession #NP_037362) fused with a Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.