

Product Datasheet

Recombinant human VEGF-C/Flt4-L/VRP Protein, His Tag (orb1551120)

Description	Recombinant human VEGF-C/Flt4-L/VRP Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	VEGFC
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1367, Recombinant human VEGF-C/Flt4-L/VRP Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Thr103-Arg227) of human VEGF-C/Flt4-L/VRP (Accession #NP_005420.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.