

## Product Datasheet

# Recombinant human Semaphorin 4A/SEMA4A Protein, Fc & His Tag (orb1551124)

<b>Description</b>	Recombinant human Semaphorin 4A/SEMA4A Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	Semaphorin 4A
<b>Note</b>	For research use only
<b>Purity</b>	> 92% by SDS-PAGE.
<b>Source</b>	KMP1363, Recombinant human Semaphorin 4A/SEMA4A Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gly32-His Tag683) of human Semaphorin 4A/SEMA4A (Accession #NP_071762.2) fused with a Fc, 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.