

Product Datasheet

Recombinant human BAFFR/TNFRSF13C Protein, Fc & His Tag (orb1551127)

Description	Recombinant human BAFFR/TNFRSF13C Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BAFFR
Note	For research use only
Purity	> 87% by SDS-PAGE.
Source	KMP1360, Recombinant human BAFFR/TNFRSF13C Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser7-Ala71) of human BAFFR/TNFRSF13C (Accession #NP_443177.1) fused with a Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.