

Product Datasheet

Recombinant human FGF acidic/FGF1 Protein, His Tag (orb1551132)

Description	Recombinant human FGF acidic/FGF1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	FGF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1355, Recombinant human FGF acidic/FGF1 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Phe16-Asp155) of human FGF acidic/FGF1 (Accession #NP_000791.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.