

Product Datasheet

Recombinant human Noggin Protein, His Tag (orb1551133)

Description	Recombinant human Noggin Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Noggin
Note	For research use only
Purity	> 85% by SDS-PAGE.
Source	KMP1354, Recombinant human Noggin Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln28-Cys232) of human Noggin (Accession #NP_005441.1) fused with a 6×His Tag at the C- terminus.
Expiration Date	6 months from date of receipt.