

## Product Datasheet

# Recombinant human Fc gamma RI/CD64 Protein, His Tag (orb1551134)

<b>Description</b>	Recombinant human Fc gamma RI/CD64 Protein, His Tag
<b>Species/Host</b>	Chinese Hamster
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	RI
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1353, Recombinant human Fc gamma RI/CD64 Protein is produced by CHO Cells expression system. The target protein is expressed with sequence (Met1-Pro288) of human Fc gamma RI/CD64 (Accession #NP_000557.1) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.