

Product Datasheet

Recombinant Mouse M-CSF/CSF-1 Protein, His Tag (orb1551136)

Description	Recombinant Mouse M-CSF/CSF-1 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	MCSF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1351, Recombinant Mouse M-CSF/CSF-1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Glu262) of mouse M-CSF/CSF-1 (Accession #NP_031804.3.) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.