

Product Datasheet

Recombinant human CEACAM-5/CD66e Protein, His Tag (orb1551138)

Description	Recombinant human CEACAM-5/CD66e Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CEACAM5
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1349, Recombinant human CEACAM-5/CD66e Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Lys35-Ala685) of human CEACAM-5/CD66e (Accession #NP_004354.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.