

Product Datasheet

Recombinant Mouse PD-L1/B7-H1 Protein, His Tag (orb1551140)

Description	Recombinant Mouse PD-L1/B7-H1 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	PDL1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1347, Recombinant Mouse PD-L1/B7-H1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Thr238) of mouse PD-L1/B7-H1 (Accession #NP_068693) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.