

Product Datasheet

Recombinant human Nectin-1/PVRL1/CD111 Protein, His Tag (orb1551142)

Description	Recombinant human Nectin-1/PVRL1/CD111 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Nectin1
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1345, Recombinant human Nectin-1/PVRL1/CD111 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Thr334) of human Nectin-1/PVRL1/CD111 (Accession #NP_002846.3) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.