

Product Datasheet

Recombinant human FGF-19 Protein, His Tag (orb1551143)

Description	Recombinant human FGF-19 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	FGF19
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1344, Recombinant human FGF-19 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Phe27-Lys216) of human FGF-19 (Accession #NP_005108.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.