

Product Datasheet

Recombinant Mouse GM-CSF/CSF2 Protein, His Tag (orb1551148)

Description	Recombinant Mouse GM-CSF/CSF2 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	GMCSF
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1339, Recombinant Mouse GM-CSF/CSF2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala18-Lys141) of mouse GM-CSF/CSF2 (Accession #NP_034099.2) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.