

Product Datasheet

Recombinant human Siglec-5/CD170 Protein, Fc & His Tag (orb1551150)

Description	Recombinant human Siglec-5/CD170 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Siglec5
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1337, Recombinant human Siglec-5/CD170 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Glu17-Thr434) of human Siglec-5/CD170 (Accession #NP_003821.1) fused with a Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.