

Product Datasheet

Recombinant human LYPD3/C4.4A Protein, His Tag (orb1551158)

| | |
|------------------------|---|
| Description | Recombinant human LYPD3/C4.4A Protein, His Tag |
| Species/Host | Human |
| Reactivity | Human |
| Endotoxins | < 0.1 EU/ug by LAL method. |
| Target | LYPD3 |
| Note | For research use only |
| Purity | > 97% by SDS-PAGE. |
| Source | KMP1329, Recombinant human LYPD3/C4.4A Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Leu31-His Tag286) of human LYPD3/C4.4A (Accession #NP_055215.2.) fused with a 6×His Tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |