

Product Datasheet

Recombinant human ADAM8/CD156a Protein, His Tag (orb1551163)

Description	Recombinant human ADAM8/CD156a Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	ADAM8
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1324, Recombinant human ADAM8/CD156a Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Pro497) of human ADAM8/CD156a (Accession #NP_001100.3) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.