

Product Datasheet

Recombinant Human/Mouse/Rat Irisin Protein, His Tag (orb1551169)

Description	Recombinant Human/Mouse/Rat Irisin Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Irisin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1318, Recombinant Human/Mouse/Rat Irisin Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Asp32-Glu143) of human Irisin (Accession #NP_715637.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.