

Product Datasheet

Recombinant human B7-H2/ICOSLG Protein, His Tag (orb1551173)

Description	Recombinant human B7-H2/ICOSLG Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	B7H2
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1314, Recombinant human B7-H2/ICOSLG Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Asp19-Ser258) of human B7-H2/ICOSLG (Accession #NP_056074.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.