

Product Datasheet

Recombinant human Cystatin SA/CST2 Protein, His Tag (orb1551174)

Description	Recombinant human Cystatin SA/CST2 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Cystatin SA
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1313, Recombinant human Cystatin SA/CST2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Trp21-Ala141) of human Cystatin SA/CST2 (Accession #NP_001313.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.