

## Product Datasheet

# Recombinant human LTBR/TNFRSF3 Protein, Fc & His Tag (orb1551176)

<b>Description</b>	Recombinant human LTBR/TNFRSF3 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	LTBR
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1311, Recombinant human LTBR/TNFRSF3 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln31-Met227) of human LTBR/TNFRSF3 (Accession #NP_002333.1) fused with an Fc, 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.