

## Product Datasheet

# Recombinant human EPCR/CD201/PROCR Protein, His Tag (orb1551182)

<b>Description</b>	Recombinant human EPCR/CD201/PROCR Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	EPCR
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1305, Recombinant human EPCR/CD201/PROCR Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser18-Ser210) of human EPCR/CD201/PROCR (Accession #NP_006395.2) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.