

Product Datasheet

Recombinant human Annexin A5/Annexin V/ANXA5 Protein, His Tag (orb1551184)

Description	Recombinant human Annexin A5/Annexin V/ANXA5 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Annexin A5
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1303, Recombinant human Annexin A5/Annexin V/ANXA5 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Asp320) of human Annexin A5/Annexin V/ANXA5 (Accession #NP_001145.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.