

Product Datasheet

Recombinant human SOST/Sclerostin Protein, His Tag (orb1551185)

Description	Recombinant human SOST/Sclerostin Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	SOST
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1302, Recombinant human SOST/Sclerostin Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln24-Tyr213) of human SOST/Sclerostin (Accession #NP_079513.1) fused with a 8×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.