

## Product Datasheet

# Recombinant Mouse HER2/ErbB2 Protein, His Tag (orb1551187)

<b>Description</b>	Recombinant Mouse HER2/ErbB2 Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	HER2
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1300, Recombinant Mouse HER2/ErbB2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Thr23-Thr653) of mouse HER2/ErbB2 (Accession #NP_001003817.1) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.