

## Product Datasheet

# Recombinant human LILRB4/CD85k/ILT3 Protein, Fc & His Tag (orb1551188)

<b>Description</b>	Recombinant human LILRB4/CD85k/ILT3 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	LILRB4
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1299, Recombinant human LILRB4/CD85k/ILT3 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln22-Glu259) of human LILRB4/CD85k/ILT3 (Accession #NP_001265356.2) fused with an Fc, 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.