

Product Datasheet

Recombinant Mouse HER2/ErbB2 Protein, Fc & His Tag (orb1551193)

Description	Recombinant Mouse HER2/ErbB2 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	HER2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1294, Recombinant Mouse HER2/ErbB2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Thr23-Thr653) of mouse HER2/ErbB2 (Accession #NP_001003817.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.