

## Product Datasheet

# Recombinant human FGF R1/CD331 Protein, Fc & His Tag (orb1551196)

<b>Description</b>	Recombinant human FGF R1/CD331 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	FGF R1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1291, Recombinant human FGF R1/CD331 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Arg22-Glu376) of human FGF R1/CD331 (Accession #NP_075598.2) fused with an Fc, 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.