

Product Datasheet

Recombinant Mouse IL-17F Protein, His Tag (orb1551197)

| | |
|------------------------|--|
| Description | Recombinant Mouse IL-17F Protein, His Tag |
| Species/Host | Human |
| Reactivity | Mouse |
| Endotoxins | < 0.1 EU/ug by LAL method. |
| Target | IL17F |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE. |
| Source | KMP1290, Recombinant Mouse IL-17F Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1- Ala161) of mouse IL-17F (Accession #NP_665855.2) fused with a 6×His Tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |