

Product Datasheet

Recombinant human STAT1 Protein, His Tag (orb1551200)

Description	Recombinant human STAT1 Protein, His Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	STAT1
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1287, Recombinant human STAT1 Protein is produced by Baculovirus-Infected Sf9 Cells expression system. The target protein is expressed with sequence (Met1- Val712) of human STAT1 (Accession #NP_644671.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.