

Product Datasheet

Recombinant Human SLAMF3 Protein, His Tag (orb1551204)

Description	Recombinant Human SLAMF3 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	SLAMF3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1283, Recombinant Human SLAMF3 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Lys48-Lys454) of human SLAMF3 (Accession #Q9HBM7) fused with a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.