

Product Datasheet

Recombinant Human BCMA/TNFRSF17 Protein, Avi & His Tag (orb1551206)

Description	Recombinant Human BCMA/TNFRSF17 Protein, Avi & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BCMA
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1281, Recombinant Human 4-1BB Ligand/TNFSF9 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Met1-Ala54) of human BCMA/TNFRSF17 (Accession #Q02223) fused with a 6xHis Tag, Avi tag at the C-terminus.
Expiration Date	6 months from date of receipt.