

## Product Datasheet

# Recombinant Human Decorin Protein, His Tag (orb1551210)

<b>Description</b>	Recombinant Human Decorin Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	Decorin
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1277, Recombinant Human Decorin Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Gly17-Lys359) of human Decorin (Accession #P07585) fused with a 6xHis Tag at the C- terminus.
<b>Expiration Date</b>	6 months from date of receipt.