

Product Datasheet

Recombinant Human 4-1BB Ligand/TNFSF9 Protein, Fc His Tag (orb1551227)

Description	Recombinant Human 4-1BB Ligand/TNFSF9 Protein, Fc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	41BB
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1260, Recombinant Human ErbB2/Her2 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Arg71-Glu254 Trimer) of human 4-1BB Ligand/TNFSF9 (Accession #P41273) fused with an Fc tag at the N-terminus.
Expiration Date	6 months from date of receipt.