

Product Datasheet

Recombinant Human TMPO Protein, His Tag (orb1551228)

Description	Recombinant Human TMPO Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TMPO
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1259, Recombinant Human TMPO Protein is produced by E.coli expression system. The target protein is expressed with sequence (Pro2-Glu187) of human TMPO (Accession #P42167) fused with a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.