

Product Datasheet

Recombinant Human Stathmin Protein, His Tag (orb1551234)

Description	Recombinant Human Stathmin Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Stathmin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1253, Recombinant Human Stathmin Protein is produced by E.coli expression system. The target protein is expressed with sequence (Ala2-Asp149) of human Stathmin (Accession #P16949) fused with a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.