

## Product Datasheet

# Recombinant Human VAMPB Protein, His Tag (orb1551235)

<b>Description</b>	Recombinant Human VAMPB Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	VAMPB
<b>Note</b>	For research use only
<b>Purity</b>	> 90% by SDS-PAGE.
<b>Source</b>	KMP1252, Recombinant Human VAMPB Protein is produced by E.coli expression system. The target protein is expressed with sequence (Ala2-Pro132) of human VAMPB (Accession #O95292) fused with a 6xHis Tag at the C- terminus.
<b>Expiration Date</b>	6 months from date of receipt.